



Manufacturer of **custom designed**

Filling & Packaging machines for liquids and semi-liquids.

Pharmaceutical

Cosmetic

Chemistry

Customized to suit your production requirements.

Manufacturer of machinery and equipment for the pharmaceutical industry since 1990, our products are designed to be integrated into existing lines or to operate independently.

Our modules are designed to be assembled together to create a complete, turnkey packaging line.



COMPLETE LIGNE - INTEGRATED OR STAND-ALONE MACHINES - REVAMPING

Filling Module

Capping Module

Labeling-Marking Module

Control Module

Integrated or stand-alone Module

Compatible with Pharma industry (Isolators, RABS, single use, etc.)

VIALS

SYRINGES

BOTTLES

SPRAY

NEST

CARTRIDGES

OTHERS

FEEDING - FILLING - CAPPING - IDENTIFICATION - CONTROL - CONVEYING

Our expertise enables us to master a wide range of skills.

Feeding:	Nest			Bulk
Filling:	Volumetric	Weight	Flowmeter	Peristaltic
Capping:	Screwing	Crimping		Pushing
Identification:	Labeling	Inkjet Printing		Laser Printing
Control:	Vision			Weight
Conveying:	Mechanical			Robotized





Discover our other products and services
on www.spc.fr

Machine manufacturer

Equipment manufacturer



Manufacturer of dosing equipment

Volumetric Filling pumps (CIP/SIP, ...)
Filling, Nitrogen and Washing needles
Other tooling

Worldwide representation through our
sales **distribution network**



They trust us



sanofi



GSK

SPC GROUP
104 Rue de Lisbonne
69140 - Rillieux-la-Pape

+33 2 76 57 01 70
sales@spc.fr