




Manufacturing of high precision tooling for liquid filling



SPC Manufacturing,
a subsidiary company of



YOUR PRECISION PARTNER
FOR LIQUID FILLING



FILLING PUMPS

PHARMACEUTICALS
R&D
CLINICAL TRIALS
BIOPHARMA
ADVANCED THERAPIES

Seal-less volumetric piston pumps are designed for packaging liquids. They can be used on all filling machines brands.

FEATURES

- Custom designed
- Monobloc stainless steel or ceramic design, to meet pharmaceutical requirements
- 316L Stainless steel
- ZrO₂ + Y₂O₃ Yttria-stabilized zirconia (white)
- ZrO₂ + MgO - stabilized with Magnesia (white or yellow)
- Al₂O₃ Alumina (>99.5%)
- Hard chrome and PVD CrN coating (providing a high level of mechanical resistance)
- CIP/SIP version available
- Watertightness-test prior to shipment (Flotest[®])
- CGMP compliant
- Decontamination and passivation
- Traceability documentation supplied
- Maintenance and repairing of existing devices



104 rue de Lisbonne, Parc d'Activités de SERMENAZ
69 140 RILLIEUX LA PAPE - FRANCE

T : +33 (0)4 78 88 79 56 • F : +33 (0)4 78 88 75 72
spc@spc.fr • <http://www.spc.fr>

